Spengler

DEPUIS 1907

ES-60®

TENSIOMÈTRE ÉLECTRONIQUE PROFESSIONNEL
AUTOMATIC BLOOD PRESSURE MONITOR
TENSIOMETRO ELECTRÓNICO PROFESIONAL
SIGMOMANOMETRO ELETTRONICO PROFESSIONALE
AUTOMATISCHE PROFI-BLUTDRUCKMESSGERÄT
AUTOMATISCHE PROFESSIONELE BLOEDDrukKETER
جهاز إلكتروني مهني لقياس ضغط الدم

Manuel d'utilisation - Instruction manual - Manual de uso
Manuale d'uso - Bedienungsanleitung - Gebruiksaanwijzing
دليل الاستعمال

WWW.SPENGLER.FR

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A HISTORIC OCCASION

In 1907, two university professors, Henri Vaquez and Charles Laubry, invented the world's first aneroid blood pressure monitor in Paris. Together with Mr Emile Spengler, they created Spengler to develop this invention.

For nearly a century, Spengler has manufactured medical instruments that enable people who use them to better understand their bodies and to live well. Spengler is today a worldwide leader in health monitoring and medical diagnostic equipment.

Spengler's development is based on a vision that includes access to better healthcare for all and instruments of superlative quality. The manufacturing processes of Spengler, which can be compared to those of watch making, are dictated by their well-known zero defect policy that has made Spengler a medical reference during the last 100 years.

Type BF applied part

Read the instructions carefully before using this device

The Spengler **ES-60®** blood pressure monitor is a fully automatic, digital blood pressure measuring device with integrated **MAM** technology for use on the upper arm. This monitor enables fast and reliable measurement of the systolic and diastolic blood pressure as well as the pulse by use of the oscillometric method.

This device is easy to use and is proven in clinical studies to provide excellent accuracy. The measurement algorithm used in this device has been validated with an A/A rating according to the strict protocol of the **BHS** (British Hypertension Society) at the St-Thomas Hospital in London*. The large display makes it easy to visualize operational status and cuff pressure during a measurement.

Please read through this instruction manual carefully before use and then keep it in a safe place.

ELECTROMAGNETIC INTERFERENCE

This device contains sensitive electronic components (micro-processor). Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens).

Electrical interference can lead to temporary impairment of the measuring accuracy.


1. MAM TECHNOLOGY


«MAM technology» is a patented technology designed for optimum reliability in measurement of blood pressure. Advanced measurement accuracy is achieved by the automatic analysis of three successive measurements. The three readings are first

* BTO-A: AN ACCUATE AUTOMATED DEVICE FOR USE IN ADULTS AND PREGNANCY, Reinders A, Cuckson AC, Lee J, Shennan AH, Maternal and Fetal Research Unit, GKT School of Medicine, St.Thomas' Hospital, London, UK

compared to one another and, depending on the relation between the three readings; the data is adjusted to give a comprehensive result rather than a simple average.

2. COMPONENTS OF THE ES-60®

3. SPENGLER PROFESSIONAL CUFFS

The cuffs are available in the different sizes (S, M, L), made of nylon. Nylon cuff can be machine-washed and spray disinfected.

4. USING THE ES-60®


4.1. Inserting the batteries

Insert and plug in battery pack immediately after unpacking the unit. The battery compartment is located under the device.

If a battery warning  appears in the display, the battery pack is discharged and must be recharged; and if a battery warning  appears, the battery pack doesn't work till being charged again.

4.2. Setting date and time

1. After the battery pack is fitted, the year number flashes in the display. You can set the year by pressing the **M**-button (or withdraw by pressing **O/I** button directly). To confirm and then set the month, press the **Time/MAP** button.

20 09

2. You can now set the month using the “M” button. To confirm and then set the day, press the **Time/MAP** button.

3. Please follow the instructions above to set the day, hour and minutes.
4. Once you have set the last minute and pressed the time button, the date and time are set and the time is displayed.

5. If you want to change the date and time, press **O/I** button and then hold the **Time/MAP** button down for approx. 3 seconds once “0” reset is done, until the year number starts to flash. Now you can enter the new values as described above.

4.3. Using the AC adapter/battery charger

This blood-pressure monitor uses an AC adapter (output 7.5 V DC / 2.0A with DIN plug). It is designed specifically for this instrument to be used as an AC/DC adapter and as a battery charger.

Please use an approved Spengler battery pack and AC/DC adapter. The use of an AC/DC adapter with the incorrect voltage or amperage will permanently damage the device.

5. CARRYING OUT A MEASUREMENT

Notes to read before using this device:

- Only use clinically approved Spengler cuffs.
- Average Mode MAM: clinical studies demonstrate taking multiple bloodpressure readings and calculating an average is more likely to determine your true blood pressure. This unit allows you to switch the unit to a special Average Mode setting (**MAM**) that automatically takes multiple readings.

If you would like to take an Average Mode measurement, please slide the switch to the right toward no. 3. A «**MAM**» symbol in the display indicates that the unit is set to the Average Mode.

When the unit is set to Average Mode, 3 separate measurements will take place in succession. There will be a 15 seconds resting time between each measurement. A count-down indicates the remaining time and a beep will sound 5 seconds before the 2nd and 3rd readings will begin. If the single data of each cycle differ too much from each other, a fourth measurement is performed before the result is displayed. In rare cases the blood pressure is so unstable that even after four measurements there is not enough data to give an accurate reading. In this case «**ERR 6**» is shown and no result can be given.


If one of measurement causes an error message it is repeated.

- Single mode:

If you would like to take single measurement, please slide the switch toward no. 1.

Taking a measurement

- Sit patient down, position unit on the table in front of you and plug in the correct sized cuff on the side of the unit.
- Fit the cuff on the arm so that the tube is facing downwards in the bend of the elbow. Gently close the cuff.

- Press **O/I** button

Select the measuring mode (**MAM** or single).

In **MAM** mode, the blood pressure will be taken three times consecutively at intervals of 15 seconds and an intelligent average will be shown.

- Press the «**START**» button (if at anytime you need to discontinue a reading, simply press the «**O/I**» button).

- After the measurement is completed, the systolic and diastolic blood pressures and the pulse rate appear in the display.

- By pressing the «**Time/MAP**» button, it is possible to display the **MAP**: Mean Arterial Pressure.

- The measurement readings remain on the display until you switch off the instrument. If no button is pressed for a period of 5 minutes the instrument switches itself off in order to preserve the batteries.

- To get the latest measurements, press on the «**M**» button.

How not to store the reading

Press the **O/I** button while the reading is being displayed. Keep the button pressed until **CL** is flashing and then release it. Confirm by pressing the **M**-button.

Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the «**O/I**» button can be pressed at any time. The pressure in the cuff is immediately reduced.

Recalibration («**TEST**» mode) in case of absurd results of measurement, press on the button «**TEST**» to recalibrate. An aural signal announces in the end of the calibration.

6. CLAIMS OF USE

6.1. Planned uses

- Measure systolic and diastolic blood pressure, pulse on arm
- Use on adults only
- Use at home and in support of treatment.

6.2. Restrictions of use

- Do not use on children or infants.
- Not suitable for measurement on an unscheduled body part, such as wrist or leg.
- Not suitable if the cuff is not placed at the correct height on the arm.
- Not suitable if the cuff used is not the right size.
- Sit in a quiet environment for at least 5 minutes taking the measurement.
- In no case should you change the doses of medication prescribed by your doctor.
- The measurement of blood pressure can not be done outside the expected environmental ranges: 10-40 ° C and relative humidity between 15-90%.
- Not suitable for use in a vibration emitting environment, electromagnetic field, direct sunlight, etc.

Thus, the ES-60® Electronic Blood Pressure Monitor is not suitable for use during patient transport and is not compatible with electrosurgery.

7. DATA MEMORY

At the end of a measurement, this instrument automatically stores each result, including date and time.

7.1. Viewing the stored values

This blood pressure monitor automatically stores the last 60 measurements values. By pressing the **M** button after **O/I** initiated and “0” reset is done, the last measurement (MR60) as well as the previous 59 measurements (MR59, MR58, ...MR1) can be displayed one after the other. Data obtained in Average Mode are indicated by a **MAM** symbol. Pressing the **M**-button repeatedly enables you to toggle between one stored value and another.

Note: press the «Time/MAP» button during memory display to check the MAP: Mean Arterial Pressure.

7.2. Memory full

Pay attention that the maximum memory capacity of 60 is not exceeded. When the memory is full, the old values are automatically overwritten with new ones. Values should be evaluated before the memory capacity is reached - otherwise data will be lost.

7.3. Clear all values

If you are sure that you want to permanently remove all stored values, hold down the **M**-button after **O/I** initiated and “0” reset is done, until **CL** appears and then release the button. To permanently clear the memory, press the **M**-button while **CL** is flashing. Individual values cannot be cleared.

8. ERROR MESSAGES/MALFUNCTIONS

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed.

Error No.	Possible causes
ERR 1	<ul style="list-style-type: none"> The systolic pressure was determined but afterwards the cuff pressure fell below 20 mmHg. The tube may have become unplugged after the systolic blood pressure was measured. No pulse could be detected
ERR 2	<ul style="list-style-type: none"> Unnatural pressure impulses impair the measurement result. Possible cause: the arm was moved during the measurement (artefact).
ERR 3	<ul style="list-style-type: none"> Inflation of the cuff takes too long. The cuff is not correctly placed or the cuff connection is not sealed.
ERR 5	<ul style="list-style-type: none"> The measured readings indicated an unacceptable difference between systolic and diastolic pressures. Take other reading following directions carefully.
ERR 6	<ul style="list-style-type: none"> Single data differ too much during an Average Mode reading even after 4 cycles. No average result can be displayed.
ERR 7	<ul style="list-style-type: none"> When «ERR 7» appears, the device is put automatically in «TEST» mode

If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

	Malfunction	Remedy
EN	The display remains empty when the instrument is switched on. Batteries have been installed.	1. Check that the battery pack is charged. 2. If the display is erratic or unusual, remove the batteries and re-install new batteries. 3. Battery pack might be damaged and needs to be replaced.
ES	The pressure does not rise although the pump is running.	1. Check the connection of the cuff tube and connect properly if necessary.
IT	The instrument frequently fails to measure the blood pressure values, or the values measured are too low (too high).	1. Ensure the cuff fits correctly. 2. Make sure that the cuff is not too tight. Ensure tight clothing such as a a rolled-up sleeve isn't exerting pressure on the arm above the measuring position. Take off articles of clothing if necessary.
DE		

Further information

For licensing, this device has been subjected to strict clinical tests and supervised by experienced Cardiac Specialists. This device is manufactured to meet European standards under the supervision of the Technical Monitoring Association. (RWTÜV-Essen).

If you have any questions regarding the use of this blood pressure monitor, please ask your dealer for the Spengler Service representative in your country. The Spengler Service Team will be happy to help you.

9. CARE AND MAINTENANCE, RECALIBRATION

- a)** Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- b)** The cuff contains a sensitive air-tight bubble. Handle this carefully.
- c)** Clean the device with a soft, dry cloth. Do not use petrol, thinners or similar solvent.
The cuff must not be immersed.
- d)** Do not drop the instrument or treat it roughly in any way. Avoid strong vibrations.
- e)** Never attempt to repair the instrument yourself!
Any unauthorized opening of the instrument invalidates all guarantee claims!

Periodical re-calibration

Sensitive measuring devices must be checked for accuracy from time to time. We recommend a periodical inspection of your unit by an authorized dealer every 2 years. Your Spengler dealer would be pleased to provide more information regarding calibration according to the EN 1060-3 / 9.2 standards.

Disposal

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

10. GUARANTEE

This professional blood-pressure monitor is guaranteed for 3 years from date of purchase. This guarantee includes the instrument but not the cuff or battery pack, which are subjected to everyday wear and tear. The guarantee does not apply to damage caused by improper handling, damage from leaking batteries, accidents, not following the operating instructions or alterations made to the instrument by third parties.

Please read the following conditions of guarantee:

Duration: 3 years from the date of purchase. The guarantee should be registered by our After-Sales Service (guarantee card to be completed and returned to us within 15 days from purchase).

FR

EN

ES

IT

DE

NL

AR

Guarantee (material and workmanship): free replacement of approved faulty parts by our After-Sales Service. The replacement or repair during the guarantee period do not extend its duration.

Consumer's statutory rights: these guarantee condition statements do not affect the consumer's statutory rights.

Guarantee exclusions: rubber and P.V.C. parts / damages due to transport / fair wear and tear of the product / damages and breakdowns due to shock, misuse or neglect / any product disassembled, repaired or modified by the user or unauthorized persons.

Guarantee service:

Warning: we advise to use your electronic BPM ES-60® for spot-checking only. We do not recommend using ES-60® in ICU or in conditions where the patient needs to be monitored.

For any intervention during the guarantee period, apply to one of the authorized Spengler After-Sales Dealers. Please contact our Customer Department for the list of Spengler representatives in your country. The guarantee applies when the product is returned undisassembled, carefully packed and postage paid, to the authorized Spengler After-Sales Dealer.

Please, send the product together with the shipment notice, a brief explanatory note and the guarantee card, if valid.

11. REFERENCE TO STANDARDS

Device standard:

Device corresponds to the requirements of the European standard for non-invasive blood-pressure monitor EN1060-1, EN1060-3, EN1060-4, DIN 58130, NIBP - clinical investigation, US standard: ANSI / AAMI SP10, NIBP - requirements

Electromagnetic compatibility:

Device fulfills the stipulations of the European standard EN 60601-1/-1-2

Clinical testing:

Clinical performance testing was carried out in Germany according to the DIN 58130 / 1997 procedure N6 (sequential).

The stipulations of the EU-Guidelines 93/42/EEC for Medical Products Class IIa have been fulfilled.

British Hypertension Society

12. TECHNICAL SPECIFICATIONS

Weight	770 g (with batteries)
Size	110 x 177 x 175 mm
Storage temperature	-20 to +55°C
Humidity	15 to 90% relative humidity maximum
Operation temperature	10 to 40°C
Display	LCD-Display (Liquid Crystal Display)
Measuring method	oscillometric
Pressure sensor	capacitive
Measuring range	
SYS/DIA	30 to 280 mmHg
Pulse	40 to 200 beats per minute
Cuff pressure display range	0-299 mmHg
Memory	Auto stores the last 60 measurements
Measuring resolution	1 mmHg
Accuracy	Pressure within \pm 3 mmHg Pulse \pm 5 % of the reading
Power source	DC 7.5V/2.0A AC adapter (DC 7.5V/2.0A) or battery pack (4.8 V DC)
Reference to standards	IEC 60601-1 / IEC 60601-1-2 EN 1060-1 / -3/-4, NIBP - requirements DIN 58130, NIBP-clinical investigations ANSI/AAMI SP10, NIBP - requirements

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SPENGLER SAS