MANUEL UTILISATEUR
USER'S MANUAL
MANUAL
دليل الملا

DEPUIS 1907

MAXI +3

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USING THE MAXI+3 LARGE DIAL ANEROID SPHYGMOMANOMETER

Steps

- Leave the patient to rest for 5 minutes
- Place the patient in a comfortable position, with legs uncrossed, and the back and arm elevated, ideally without moving or speaking
- Check that the marker index for arterial index is within the range shown on the Maxi+3 cuff
- Keep the cuff at heart height during measurement
- Keep the pneumatic circuit closed (check the decompression valve)
- Inflate the cuff for an initial SBP estimation (when the blood circulates again, it causes a noticeable beating sound on the stethoscope) until radial pulse disappears (so-called Korotkoff sounds)
- Inflate again to 50mmHg above the estimated SBP
- Deflate by opening the valve at a rate of 2mmHg per second and simultaneously read the blood pressure on the manometer
- Read the SBP for phase I of Korotkoff sounds
- Read the DBP for phase V of Korotkoff sounds (complete disappearance of sounds noted at the stethoscope)
- Measure the heart rate
- Continue decompression until complete purge of the cuff


⚠ Warning: A blood flow stoppage carries risks.

Recommandation:

- Use K5* for adult auscultation
- Use K5* for auscultation of pregnant patients, so long as the sounds are noticeable when the cuff is deflated, in which case the K4* can be used
- Use K4* for auscultation of children between 3 and 12 years.

*See the table describing the phases of Korotkoff sounds.

KOROTKOFF SOUNDS

ENVIRONMENTAL SPECIFICATIONS:

Operating temperature: +10°C to +40°C
Operating relative humidity: 20% to 85%
Storage temperature: -20°C to +70°C
Relative humidity at storage: 20% to 85%

MAINTENANCE AND MAINTENANCE VERIFICATION

Comply with current regional regulations for the disposal of manual tensiometers or accessories.

METROLOGICAL INSPECTION BY AN APPROVED BODY

Check accuracy and watertightness:

- At least once every 2 years (manometer)
- At least once per year (cuffs, tube, bulb)
- After each maintenance operation or repair or replacement of removable parts.

Note: In accordance with international and national metrological standards, the operator or the body responsible for maintenance must ensure traceability of benchmark manometer used for calibration.

REMOVABLE PARTS & ACCESSORIES:

Please contact your approved retailer to purchase the parts and accessories below:

- Replacement switch,
- Tube,
- Armbands (all sizes).

CLEANING AND DISINFECTION:

⚠ Warning: Never immerse the manometer in liquid or in a decontamination device. Use a neutral detergent to clean the surface of product. Should disinfection be necessary, use a cloth soaked in 75% alcohol.

⚠ Warning: Avoid placing your sphygmomanometer too close to heat, excessive cold, solvents or greases. Avoid long-term contact with sweat, as the PVC could easily stiffen (switch and tube). Keep hands dry.

TECHNICAL SPECIFICATIONS:

Maxi+3 Manometer

Accuracy: +/-3mmHg

Operating range: 0 to 300mmHg

Displayed range: 0 to 300mmHg

Indication type: radial scale

Scale grading: increments of 2 mmHg

Pressure generation: switch

Pressure reduction: air release valve

Case: ABS

Mechanism: copper-beryllium alloy

Switch: PVC

Tube: PVC.

Cuff: Nylon.

PRECAUTIONS FOR USE

Avoid impacts, particularly on the dial glass.

Do not expose to heat, sun or humidity, dust.

WARRANTY

The MAXI+3 sphygmomanometer is guaranteed against any part and labour or operating fault, in accordance with the manufacturer's specifications, under normal maintenance and usage conditions. The 2-year guarantee for the tensiometer (excluding armband and tube) starts from the date of purchase from an approved retailer. SPENGLER's obligations are limited during the guarantee period to the repair or replacement of components deemed to be defective by SPENGLER.

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