Spengler

TEMPO EASY

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EN

TEMPO EASY

Non-contact Infrared Thermometer
User's manual

IMPORTANT

Read this manual completely before using the Spengler non-contact Infrared Thermometer.
Save these instructions.

CONTENTS

EN

1. Safety precautions	p.26
2. Introduction	p.27
3. Precautions before use	p.27
4. Operating principle.....	p.28
The different methods of temperature measurement	p.28
Normal temperatures according to measurement method	p.29
Advantages of temporal artery (TA) temperature	p.30
Normal temperature according to age	p.30
Practical considerations when taking a temperature	p.30
How to take a temperature	p.31
Constraints	p.32
5. Basic instrument.....	p.33
6. Features	p.34
7. Instructions	p.34
8. Setting and function of menu	p.35
9. Technical specifications.....	p.37
10. Maintenance of the product	p.38
11. Accessories.....	p.38
12. Guidelines	p.39
13. Troubleshooting	p.45

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1.SAFETY PRECAUTIONS

- Follow the maintenance advice stipulated in this instruction manual.
- This device may be used for professional purposes or for personal home use.
- This device must only be used for the purposes described in this instruction manual.
- This device must only be used in an ambient temperature range of between 10°C and 40°C.
- This device must always be kept in a clean, dry area.
- Do not expose this thermometer to electric shocks.
- Do not expose this thermometer to extreme temperature conditions of 55°C or -20°C.
- Do not use this device in relative humidity higher than 85%.
- The protective glass over the lens is the most fragile part of the thermometer.
- Do not touch the glass of the infrared lens with your fingers.
- Clean the glass with a cotton bud lightly moistened with 95° alcohol.
- Do not expose the thermometer to sunlight or to water.
- Never drop the device.
- Should a problem occur with your device, please contact your retailer.
- Do not attempt to repair this device yourself.

THE MANUFACTURER RESERVES THE RIGHT TO
ALTER THE SPECIFICATIONS OF THE PRODUCT
WITHOUT PRIOR NOTIFICATION

2.INTRODUCTION

The TEMPO EASY Non-contact Infrared thermometer has been developed by using the latest infrared technology. This technology allows temporal artery (TA) temperature to be taken at a distance of about 3cm~5cm away from the forehead. Precise, Instantaneous and without Contact, the Tempo Easy is, up to now, the most suitable thermometer for no risk on temperature measurement. It has been demonstrated that this method of TA temperature measurement is more precise than the tympanic thermometry and better tolerated than rectal thermometry (1). However, as with other types of thermometer, it is essential to use the Tempo Easy properly in order to obtain reliable and stable results. You are therefore advised to read this instruction manual and the safety precautions carefully before use.

(1)Greenes D, Fleisher G. Accuracy of a Noninvasive Temporal Artery Thermometer for Use in Infants. Arch Pediatr Adolesc Med 2001;155:376.

3.PRECAUTIONS BEFORE USE

The TEMPO EASY is pre-set at the factory.
It is not necessary to calibrate the device when starting it up.

In order to obtain reliable and stable results, you are advised each time there is a significant change in the ambient temperature due to a change in environment, to allow the TEMPO EASY to acclimatize to this ambient temperature for 15 to 20 minutes before using it. It is important to allow 3~5 seconds interval between two measurements.

4. OPERATING PRINCIPLE

All objects, solid, liquid or gas, emit energy by radiation. The intensity of this energy depends on the temperature of the object. The TEMPO EASY infrared thermometer is therefore able to measure the temperature of a person by the energy the person emits. This measurement can be taken thanks to an external temperature probe on the device which permanently analyses and registers the ambient temperature. Therefore, as soon as the operator holds the thermometer near the body and activates the radiation sensor, the measurement is taken instantly by detection of the infrared heat generated by the arterial blood flow. Body heat can therefore be measured without any interference from the heat of the surrounding environment.

THE DIFFERENT METHODS OF TEMPERATURE MEASUREMENT

Core temperature

Core temperature is the most precise measurement and involves measuring the temperature in the pulmonary artery by means of a catheter equipped with a thermal probe which can read the temperature in situ. The same method is employed for probes measuring the oesophageal temperature. However, such invasive temperature measurement methods require specific equipment and expertise.

Rectal thermometry

Rectal temperature adjusts slowly in comparison to the evolution of the body's internal temperature. It has been demonstrated that rectal temperature remains raised long after the internal temperature of the patient has started to drop and vice versa.

Furthermore, rectal perforations have been known to occur as a result of this method and without appropriate sterilisation techniques, rectal thermometry can spread germs often found in faeces.

Oral thermometry

Oral temperature is easily influenced by recent ingestion of food or drinks and by breathing through the mouth. To measure oral temperature, the mouth must remain closed and the tongue lowered for three to four minutes which is a difficult task for young children to accomplish.

Axillary (armpit) temperature

Although it may be easy to measure axillary temperature, it has been proven that it does not provide an accurate measurement of the child's internal temperature. To take this type of temperature, the thermometer must be wedged tightly over the axillary artery. Despite the low sensitivity and relative inaccuracy of axillary temperature in detecting fever, this method is recommended by The American Academy of Pediatrics as a screening test for fever in newborns.

Tympanic thermometry

In order to obtain a precise temperature reading, good command of the measurement technique is required. The thermometer probe must be placed as close as possible to the warmest part of the external ear canal. An incorrectly placed probe could lead to a false temperature reading.

NORMAL TEMPERATURES ACCORDING TO MEASUREMENT METHOD

MEASUREMENT METHOD	NORMAL TEMP°
RECTAL	36.6°C – 38°C
ORAL	35.5°C – 37.5°C
AXILLARY	34.7°C – 37.3°C
AURICULAY	35.8°C – 38°C
TEMPORAL	35.8°C – 37.8°C

The temperature of the human body varies throughout the day. It can also be influenced by numerous external factors: age, sex, type and thickness of skin...

ADVANTAGES OF TEMPORAL ARTERY (TA) TEMPERATURE

Infrared arterial temperature can be measured using a device placed on the forehead, in the temporal artery region. It has been demonstrated that this relatively new method of measuring temperature is more precise than tympanic thermometry and better tolerated than rectal thermometry.

The TEMPO EASY thermometer has been designed to produce an instant forehead temperature reading without any contact with the temporal artery. As this artery is quite close to the surface of this skin and therefore accessible and given the blood flow is permanent and regular, it allows precise measurement of the temperature. This artery is linked to the heart by the carotid artery which is directly linked to the aorta. It forms part of the main trunk of the arterial system. The efficiency, speed and comfort of taking a temperature from this area make it ideal compared with other temperature measurements methods.

NORMAL TEMPERATURES ACCORDING TO AGE

Age	°C	°F
0-2 ans	36,4-38,0	97,5-100,4
3-10 ans	36,1-37,8	97,0-100,0
11-65 ans	35,9-37,6	96,6-99,7
> 65 ans	35,8-37,5	96,4-99,5

PRACTICAL CONSIDERATIONS WHEN TAKING A TEMPERATURE

- In order to ensure that precise and accurate temperature measurements are obtained, it is essential that each user has received adequate information on and training in the temperature measurement technique when using such a device.
- It is essential to remember that although procedures such as taking a temperature may be simple they must not be trivialised.

- Temperature should be taken in a neutral context. The patient must not have undertaken vigorous physical activity prior to taking his/her temperature and the room temperature must be moderate.
- Be aware of physiological variations in temperature which must be taken into consideration when evaluating the results: temperature increases by 0.5°C between 6am and 3pm. Women have a temperature that is higher, on average, by around 0.2°C. Their temperature also varies in accordance with their ovarian cycle. It rises by 0.5°C in the second half of the cycle and at the early stages of pregnancy.
- When sitting, temperature is lower by about 0.3°C to 0.4°C than when standing.

HOW TO TAKE A TEMPERATURE

Aim at the FOREHEAD, over the right temporal region, from a distance of about 3cm~5 cm, press the thermometer's measurement button and the temperature is instantly displayed.

The reliability of the measurement cannot be guaranteed if the temperature is measured over another part of the body (e.g.arm, torso...)


CONSTRAINTS

Please observe the following before any temperature measurement to ensure a stable and reliable result:

- Push back hair from the forehead
- Wipe away any perspiration from the forehead
- Avoid any drafts (e.g. from nasal specs, air conditioning...)
- Allow a 1 minute interval between two measurements.
- Each time there is a significant change in the ambient temperature due to a change in environment, to allow the TEMPO EASY to acclimatise to this ambient temperature for at least 15 minutes before using it.

5.BASIC INSTRUMENT

EN

Symbole	Reference
	Trade mark
	IEC 60417-5333, Type BF applied part
	IEC 60417-5031 Direct current
IP22	IEC 60529 Ingress protection
	WEEE Directive (2002/96/EC)
	Consult operating instructions

6.FEATURES

1. Special design to take the Human Body Temperature with a 3 cm~5 cm (1.2~2 in) distance from forehead.
2. Reliable and stable measurement, thanks to the advantage Infrared Detection System.
3. Audible alarm if temperature is more than 38°C (100.4°F).
4. Memorize the last 32 temperature measurements.
5. Backlits LCD digital display screen.
6. Temperature unit can be displayed in either Celsius or °C
7. Automatic power-off (30 secs) to conserve energy.
8. Longevity use (100,000 readings).
9. Practical, easy to use.

ADDITIONAL USAGE


TEMPO EASY can also be used to measure the temperature of a baby-bottle or bath, or room temperature (by using the Surface Temp function).

This function is in accordance with the Directive.

7.INSTRUCTIONS

1. Install battery.
2. For the first use or when inserting new batteries wait between 10 minutes for the warm-up of the apparatus and when inserting the new batteries.
3. Aim towards the forehead (see the diagram below for the Tempo Easy positioning), from a distance of 3 cm~5 cm, press the measuring key, the temperature is displayed in 1 second.
4. Before taking the temperature, make sure to remove hair and perspiration from the forehead.

8.SETTING AND FUNCTION OF MENU

1.Switch on the device

In the switch on state, Press the “SCAN” button, one second after the screen panel in full display, it will enter the standby mode with the sign “---°C” or “---°F”.

Then press the “SCAN” button again, you will get the measuring result in 1 second. But if there is no more operation, it will turn off in 30 seconds automatically.

2.In the switch on state, Setting the mode

a. Press “MODE” button, and the screen will display:

Body...°C

b. Press again “MODE” button and the screen will display: Surface
Temp...°C

c. Press again “MODE” button and the screen will display: Room...°C

Note: The thermometer default is set to BODY mode.

3.In the switch on state, Press the “°C/ °F” button to transfer between Celsius and Fahrenheit.

4.In the switch on state, Press “MEM” (Memory) button, which will then display the last temperature, and allows for a view of the last 32 measurements.

Important!

The surface body temperature differs from the internal body temperature. To obtain the internal temperature, always use the "BODY" mode. Please make sure to select "SURFACE TEMP" mode for an external area reading.

5. In the switch on state, press “” can open or close the buzzer.

When the screen shows “ON”, the buzzer opened.

When shows “OFF”, buzzer closed.

6. Recalibration of device via the F4 MENU

The TEMPO EASY factory-set. Depending on various skin types and thickness, there may be temperature differences.

The F4 function will recalibrate and if you see a difference in the reading ,you must recalibrate the unit.

Instructions for recalibration :

In the switch on state, Press the “MODE” and “°C/ °F” buttons at the same time for 3 seconds.

The screen will then display : “F4”

Press “MODE” to increases 0.1°C, press °C/ °F to reduce 0.1°C.


Press “MEM” button to save the setting.

In the cases of seasonal or environmental changes a verification and adjustment should be carried out.

7. CHANGING THE BATTERIES

Display: when the LCD screen displays “Battery”, the battery is used.

Operation: Open the lid and change the batteries, taking great care with the correct positioning. A mistake with this could cause damage to the apparatus and compromise the guarantee of your TEMPO LASER. Never use rechargeable batteries. Use only batteries for single usage.

9.TECHNICAL SPECIFICATIONS

1.Normal using condition

Ambient temperature: 10°C ~ 40°C (50°F ~ 104°F)

Relative humidity : ≤ 85%

2.Batteries: DC 3V (2 pcs AA batteries)

3.Unit size: 175 x 50 x 52 mm (L x W x H)

4. Unit weight (including batteries): 48g

5.Temperature display resolution: 0.1°C (0.1°F)

6.Measuring range:

In body mode: 32°C ~ 43°C (89.6°F ~ 109.4°F)

In surface temp mode: 0°C ~ 60°C (32°F ~ 140°F)

In room mode: 0°C ~ 40°C (32°F ~ 104°F)

7.Precision:

32.0°C ~ 34.9°C (89.6°F ~ 94.8°F) ±0.3°C(±0.6°F)

35.0°C ~ 42.0°C (95.0°F ~ 107.6°F) ±0.2°C(±0.4°F)

42.1°C ~ 43°C (107.8°F ~ 109.4°F) ±0.3°C(±0.6°F)

8. Consumption: ≤ 300mW
9. Measuring distance: 3cm ~ 5cm (1.2in ~ 2in)
10. Automatic power-off: 30 secs
11. Memory: 32 sets

Note : The Non-contact Infrared Thermometer TEMPO EASY can take temperature readings below 32°C or above 43°C (89.6°F to 109.4°F) but precision is not guaranteed outside of this range.

LONGEVITY OF THE PRODUCT

The TEMPO EASY was conceived for an intense and professional use, its longevity is guaranteed for 40,000 takings.

10. MAINTENANCE OF THE PRODUCT

- The protective glass over the lens is the most important and fragile part of the thermometer, please take great care of it.
- Clean the glass with cotton fabric, wet with 95° alcohol.
- Do not use other batteries than mentioned batteries, do not recharge non rechargeable batteries, do not throw in fire.
- Remove the batteries when thermometer is not used for an extended period of time.
- Do not expose the thermometer to sunlight or water.
- An impact will damage the product.

11. ACCESSORIES

- User Manual 1 pc
AA alkaline batteries 2 pcs
Carry bag 1 pc

12.GUIDELINES

This device complies with the EU Directive 93/42/EC concerning medical products, the ISO 80601-2-56 and the European Standard EN60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility.

EMC DECLARATION

Guidance and manufacturer's declaration – electro-magnetic immunity			
The TEMPO EASY is intended for use in the electro-magnetic environment specified below. The customer or the user of the TEMPO EASY should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electro-magnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Not Applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to line(s)	Not Applicable	Mains power quality should be that of a typical commercial or hospital environment.

Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	Not Applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the TEMPO EASY requires continued operation during power mains interruptions, it is recommended that the TEMPO EASY be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electro-magnetic immunity

The TEMPO EASY is intended for use in the electromagnetic environment specified below. The customer or the user of the TEMPO EASY should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
			<p>Portable and mobile RF communications equipment should be used no closer to any part of the TEMPO EASY, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p>
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Not Applicable	$d = \left[\frac{3.5}{v1} \right] \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	$d = \left[\frac{3.5}{E1} \right] \sqrt{P}$ <p>80MHz to 800MHz</p> $d = \left[\frac{7}{E1} \right] \sqrt{P}$ <p>800MHz to 2.5 GHz</p>

		<p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol:</p> <p></p>
--	--	---

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the TEMPO EASY is used exceeds the applicable RF compliance level above, the Medical TEMPO EASY should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the TEMPO EASY.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Medical TEMPO EASY

The TEMPO EASY is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Medical TEMPO EASY can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the TEMPO EASY as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150KHz to 80MHz	80MHz to 800MHz	800MHz to 2,5GHz
	$d = \left[\frac{3.5}{V1} \right] \sqrt{P}$	$d = \left[\frac{3.5}{E1} \right] \sqrt{P}$	$d = \left[\frac{7}{E1} \right] \sqrt{P}$
0,01	/	0.12	0.23
0,1	/	0.38	0.73
1	/	1.2	2.3
10	/	3.8	7.3
100	/	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

13.TROUBLESHOOTING

If you have problems while using your thermometer, please refer to this guide to help resolve the problem. If the problem persists, please contact your local retailer.

THE SCREEN DISPLAYS TEMPERATURE HIGHER THAN 42.9°C (109.2°F):

The temperature is in Fahrenheit. Change the measurement to Celsius by pressing the Celsius/Fahrenheit button.

THE SCREEN DISPLAYS TEMPERATURE LOWER THAN 32°C (89.6°F):

To take the surface temperature, press the “Mode” button and set to the reading called “Body”, If the device is in Surface Temp Mode, the 32°C (89.6°F) temperature displayed is showing the external temperature of your body, rather than the internal.

THE SCREEN DISPLAYS THE MESSAGE**HI**

When using the TEMPO EASY Thermometer, the message "HI" can show on the screen. In this case, the temperature is above the measurement range selected, either above 42.9 °C (109.2°F) in Body Mode.

**THE SCREEN DISPLAYS THE MESSAGE****LO**

When using the TEMPO EASY Thermometer, the message "LO" can show on the screen. In this case, the temperature analyzed is under the measuring range selected, either less than 32°C (89.6°F) in Body Mode.

This message displays for various reasons. Please find below a list of the main issues :

Reasons for LO message display	Advice
Temperature reading hampered by hair or perspiration.	Make sure there is no obstruction or dampness prior to taking temperature.
Temperature hampered by an air draft or dramatic change in ambient temperature.	Make sure there is no air blowing in the area of use; this could affect the infrared reading.
Temperature readings are too close together, and the thermometer did not have time to reboot.	Pause for 3~5 seconds minimum between readings; a 15 seconds pause is recommended.
The measuring distance is too far.	Take measurements at the recommended distance (app. 3~5 cm; 1.2in~2.0in).

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