Spengler

OXYGO®

DB12

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Please be sure to read the manufacturer's
User Manual first and follow the instructions carefully.

Contents

CHAPTER 1	Safety Information
	1.1 Warnings 1.2 Cautions
CHAPTER 2	Introduction of Oxygo® (DB12) Pulse Oximeter
	2.1 Product Introduction 2.2 Package Contents
	2.3 Oxygo® (DB12) Pulse Oximeter Introduction
CHAPTER 3	Oxygo® (DB12) Pulse Oximeter Pre-Operational Procedures
	3.1 Battery Installation and Replacement
CHAPTER 4	Oxygo® (DB12) Pulse Oximeter Operation
	4.1 Oxygo® (DB12) Pulse Oximeter Measurement
CHAPTER 5	Oxygo® (DB12) Pulse Oximeter Cleaning and Maintenance
CHAPTER 6	Trouble shooting
CHAPTER 7	Oxygo® (DB12) Pulse Oximeter Specifications
	7.1 Technical Specifications
CHAPTER 8	Manufacturer's Declaration
CHAPTER 9	Symbol Definition
CHAPTER 10	Customer Service Information

Chapter 1 Safety Information

1.1 Warning

- Explosion Hazard : Please do not use this device in the presence of flammable gas.
- Please do not use this device in an MRI environment.
- The sole objective for the use of this device is to serve as an aid in clinical diagnosis. It must be used in combination with clinical manifestations and symptoms for the diagnosis.
- Should the subject's measurement site experiences trauma or other medical conditions resulting in incorrect result, the user shall consult the physician prior to using it.
- In the event of malfunction, please do not use this device.
- This device is able to accurately measure the pulse and blood oxygen saturation. When the blood flow is limited (such as using a sphygmomanometer), the pulse measurement may be affected. The device function may be influenced if the test object is moving.
- Do not modify this equipment without autorisation of the manufacturer.

1.2 Caution

- Not suitable for use in an oxygen rich environment.
- This is a device with a highly sophisticated electronic instrument and its maintenance and repair can only be performed by professionally trained personnel.
- Please do not expose this device to sterilization under high pressure or immerse it in liquid such as directly under the rain, as it may cause device malfunction or inaccurate measurement.
- Please refer to the cleaning procedure in this User Manual. Please do not use corrosive detergent or cleaning apparatus.
- This device is not intended for respiratory measurement.
- This device is designed to measure the ratio of the artery blood oxygen saturation to the functional hemoglobin. Functional disorders such as carboxyhemoglobinemia or methemoglobinemia may affect measurement accuracy.
- Indigo cyanine green or certain pigments in the blood vessels may have an impact on the blood oxygen concentration measured depending on its concentration.
- Defibrillators may influence the functioning of this device.
- This device may not be used under all kinds of circumstances. If unstable readings appear during measurement, please cease to use it.
- Do not swallow small component such as battery cover or battery.

- This device possesses some resistance to shaking. However, under certain circumstances, the movement may be accepted as normal pulses. Therefore, moving the patient should be avoided as much as possible.
- If the battery has not been used for a long period of time, some chemical substances may be released. If the device is to be left unused for more than 3 months, the battery should be removed.
- After the battery is installed properly, turn on the device. The device will detect the pulse signals after the finger is placed inside the device. The standby mode will be displayed within 30 seconds if no finger is placed inside the device. Please do not use it and contact the dealer for instruction if there is no display.
- The disposal and recycling of the device and its components should be in compliance with local environmental regulations.
- Under conditions of normal battery power supply, there is no limitation on the longest time of use for the **Oxygo® (DB12)** Pulse Oximeter at one single location.
- The silicon used in the **Oxygo® (DB12)** Pulse Oximeter to be in contact with the finger is a bio-compatible material and not harmful to the human body. You can therefore use it without harm.
- Keep the device out of the reach of Children, Infant, or pets.
- Some materials may cause potential allergic reaction, please stop using and find medical attention if there is any discomfort.
- Prolonged contact with skin may cause irritation, please stop using and find medical attention if there is any discomfort.
- Lint, dust, light might impact measurement.
- Some equipment might impact measurement.
Example 1 Heat from a fireplace or radiant heater might interfere the equipment.
Example 2 moisture from a nebuliser or steam kettle might interfere the equipment.
- No sterilization device

Chapter 2 Introduction of Oxygo® (DB12) Pulse Oximeter

2.1 Product Introduction

Intended Use :

The **Oxygo® (DB12)** is intended to measure blood oxygen saturation and pulse rate of adults in a non-invasive way. It is designed to be used in clinics, hospitals and home care environment.

Measurement Principle :

The measurement of SpO₂ (blood oxygen saturation concentration) is made by determining the absorption of infrared and far-infrared passing through the tissues.

Such absorptions may vary with the blood pulsation in the capillary bed allowing the blood oxygen saturation and pulse alteration to be determined.

Intended User Profile:

Patient population: Adult. Weight 47-82kg / Education: above junior high school. User group: * Users: Patient is an intended operator. *Operator: nursing staff or caregivers.

Knowledge base: Basic knowledge of human body, understand hygiene. Part of the body contacted with the device is finger.

2.2 Package Contents

The **Oxygo® (DB12)** Pulse Oximeter package includes the following :


1. **Oxygo® (DB12)** Pulse Oximeter x 1
2. User Manual
3. LR03 Alkaline Batteries x 1 (This device use internal power)
4. Pouch x 1

Please verify all of the above items are in the package. If any one item in this list is missing or damaged, please contact the distributor.

The standard configuration does not include the accessories of the system. For the accessories and the detailed information about the **Oxygo® (DB12)** Pulse Oximeter, please contact a distributor.

2.3 Oxygo® (DB12) Pulse Oximeter Instruction

1. Device Appearance and Display

2. **Oxygo® (DB12)** Pulse Oximeter Display Warning

When test data exceeds the factory default setting, the reading will flash.

Blood Oxygen Warning : When the test result falls below 90, the reading of SpO2 will flash.

Pulse Rate Warning : When the test result falls below 50 or goes above 160, the PR reading will flash.


The **Oxygo® (DB12)** has no audible alarms. It will not sound if the amount of oxygen in your blood is low or if your pulse rate is too high or too low.

Chapter 3 **Oxygo® (DB12)** Pulse Oximeter Pre-operational Procedures

3.1 Battery Installation and Replacement

This device uses one LR03 battery for its power supply. When the device displays low battery level, the battery power of **Oxygo® (DB12)** Pulse Oximeter can last only for approximately 15 minutes more. However, it may vary with different brands of batteries. Therefore, replacing the battery immediately is recommended to ensure accurate measurements.

Please follow the steps below for battery replacement.

- Step 1 Turn off power
- Step 2 Push the battery back cover down toward the end of the **Oxygo® (DB12)** Pulse Oximeter and detach the battery cover
- Step 3 Install one LR03 alkaline battery as illustrated in the figure below
- Step 4 Put the battery cover back on.

Chapter 4 Oxygo® (DB12) Pulse Oximeter Operations

4.1 Oxygo® (DB12) Pulse Oximeter Measurements

- Step 1 Press the power switch to turn on the device
- Step 2 Insert the finger
- Step 3 When the pulse signal is detected, the oxygen saturation level and pulse rate will be displayed 6~8 seconds later
- Step 4 Take the finger out of the device. The results will be displayed on the screen

Chapter 5 Oxygo® (DB12) Pulse Oximeter Cleaning and Maintenance

For the cleaning and maintenance of the **Oxygo® (DB12)** Pulse Oximeter, please follow the following recommendations:

- Prior to cleaning, turn off the **Oxygo® (DB12)** Pulse Oximeter.
 - Use a gentle detergent or medical grade alcohol (75% ethanol solution) to moisturize a swab to wipe the surface of the **Oxygo® (DB12)** Pulse Oximeter.
 - When you spot any soiling, dust, or clogging on the surface of the **Oxygo® (DB12)** Pulse Oximeter, please clean it. And it is recommended to clean once a month.
- Make sure the optical components of the device are not polluted by dust or soil.

- The same solution can be used for the cleaning and disinfection of the internal probe of the **Oxygo® (DB12)** Pulse Oximeter. Please ensure that the probe is wiped dry prior to another use. The plastic material of the probe of the **Oxygo® (DB12)** Pulse Oximeter is of medical grade. Therefore, it is not harmful or poisonous to the human skin.
- Keeping the **Oxygo® (DB12)** Pulse Oximeter dry at all times is highly recommended. Placing it in a humid environment may affect the life of the device or even cause damage to it.
- Caution: Do not spray, spill, or dump any liquid on the **Oxygo® (DB12)** Pulse Oximeter, its components, power switch or crevices.
- If the **Oxygo® (DB12)** Pulse Oximeter will not be used for a certain period of time, please take the battery out of the battery chamber.

When the low battery level is displayed, please replace the battery and comply with the local government regulations to dispose of the dead batteries.

Chapter 6 Trouble Shooting

In the event of device failure during its use, please follow the recommended solutions below.

Table 1 Recommended Breakdown Solutions

Incident	Potential Cause	Solution
Device unable to power on	No battery installed	Install battery
	Battery power depleted	Replace battery
	Battery not installed properly	Re-install battery
	Partial damage to metal dome in contact with battery	Contact distributor
No display on screen	Damaged display or connector	Contact distributor
No data displayed on device	Low perfusion	If the pulse wave is not displayed 1) Readjust finger position 2) Keep finger warm 3) Try another finger
	LED display not illuminating	Contact distributor
Automatic power off failure	Damaged electronic components	Contact distributor
No display the pulse strength	Signal inadequacy	If the pulse wave is not displayed 1) Readjust finger position 2) Keep finger warm 3) Try another finger

Chapter 7 Oxygo® (DB12) Pulse Oximeter Specifications

7.1 Technical Specifications

Model	Oxygo® (DB12)	
Display Range	SpO2	35% ~ 100%
	Pulse Rate	30 bpm ~ 250 bpm
Declared Range	SpO2	70% ~ 100%
	Pulse Rate	30 bpm ~ 200 bpm
	Perfusion Range	0.5% ~ 20%
Oxygen Saturation Accuracy	Adult	70% ~ 100% Tolerance 3% (A_{TAS})
Pulse Rate Accuracy	Adult	30 bpm ~ 200 bpm Tolerance 3 bpm (A_{PRAS})
Expected service life	Life Time	5 years
	Battery Life	About 15 hours (LR03 Alkaline Batteries)
Data update period	<3s	
Environmental Conditions	Temperature	5°C ~ 45°C (41°F ~ 113°F)
	Atmospheric Pressure	70kPa ~ 106kPa(525 mmHg ~ 795 mmHg)
	Relative Humidity	15% ~ 95% non-condensing
Shipping and Storage	Temperature	-25°C ~ +70°C (-13°F ~ 158°F)
	Atmospheric Pressure	50kPa ~ 106kPa(375 mmHg ~ 795 mmHg)
	Relative Humidity	10% ~ 95% non-condensing
Physical Characteristics	Weight	40g (battery included)
	Dimension	65mm x37mmx29mm
Electric Specifications	Electricity Type	One LR03 alkaline battery
	Voltage	1.5V DC
Compliance	Protection Category	BF Type - Application
	Operating Mode	Spot-Checking
	Device Classification Safety Standard	GB 9706.1-2007 IEC-60601-1
	Electromagnetic Compatibility	YY 0505-2005 IEC-60601-1-2, EN301489-1
	Function Standard	80601-2-61:2017
Blood Oxygen Sensor	Wavelength	660 nm, 905 nm
	Total Optical Energy Released from Probe	<15mW

The PULSE OXIMETER EQUIPMENT is calibrated to display FUNCTIONAL OXYGEN SATURATION.

Chapter 8 Manufacturer's Declaration

The **Oxygo® (DB12)** is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the monitor as recommended below, according to the maximum output of the communications equipment.

Table 5 : Recommended Separation Distances

Recommended separation distances between portable and mobile RF communications equipment and the **Oxygo® (DB12)** (IEC 60601-1-2)

Frequency of Transmitter	150 kHz ~ 80 MHz	80 MHz ~ 800 MHz	800 MHz ~ 2.5 GHz
Equation	$d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	$d = \left[\frac{7}{V_1} \right] \sqrt{P}$
Rated Maximum Output Power of Transmitter in Watts(W)	Separation Distance in Meters (m)	Separation Distance in Meters (m)	Separation Distance in Meters (m)
0.01	0.1	0.1	0.2
0.1	0.4	0.4	0.7
1	1.2	1.2	2.3
10	3.7	3.7	7.4
100	11.7	11.7	23.3

For transmitters rated at a maximum output power not listed above, the separation distance can be estimated using the equation in the corresponding column, where P is the maximum output [power rating of the transmitter in watts (W)] according to the transmitter manufacturer.

Note : At 80 MHz to 800 MHz, the separation distance for the higher frequency range applies.

Note : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Table 6 : Electronic Emissions

The **Oxygo® (DB12)** is intended for use in the electromagnetic environment specified below. The customer or user of the monitor should assure that it is used in such an environment.

Emission Test	Compliance	Rated Maximum Output Power of Transmitter in Watts(W)
RF Emissions CISPR 11	Group 1	The Oxygo® (DB12) uses RF energy only for its internal function. Therefore, the RF emissions are very low and not likely to cause interference in nearby electronic equipment. The Oxygo® (DB12) is suitable for use in establishments, including diagnostic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF Emissions CISPR 11	Class B	
Harmonic Emissions IEC 61000-3-2	N/A	
Harmonic Emissions IEC 61000-3-2	N/A	

Table 7 : Electromagnetic Immunity

The **Oxygo® (DB12)** is intended for use in the electromagnetic environment specified bellow. The customer or user of the monitor should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level		Electromagnetic Environment Guidance	
Conducted RF IEC 61000-4-6	3 Vrms ; 6 Vrms 150 kHz to 80 MHz	N/A		Portable and mobile RF communications equipment should be used no closer to any part of the EQUIPMENT or SYSTEM including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Interference may occur in the vicinity of equipment marked with the following symbol. 	
Radiated RF IEC 61000-4-3	3 V/m ; 10V/m 80 MHz – 2.7 GHz 80%	3 V/m ; 10V/m 80 MHz – 2.7 GHz 80%			
Proximity fields from RF wireless Communications equipment IEC 61000-4-3	27 V/m	385 MHz	27 V/m		385 MHz
	28 V/m	450 MHz	28 V/m		450 MHz
	9 V/m	710 MHz	9 V/m		710 MHz
		745 MHz			745 MHz
		780 MHz			780 MHz
	28 V/m	810 MHz	28 V/m		810 MHz
		870 MHz			870 MHz
		930 MHz			930 MHz
	28 V/m	1720 MHz	28 V/m		1720 MHz
		1845 MHz			1845 MHz
		1970 MHz			1970 MHz
	28 V/m	2450 MHz	28 V/m		2450 MHz
	9 V/m	5240 MHz	9 V/m		5240 MHz
		5500 MHz			5500 MHz
		5785 MHz			5785 MHz

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5 kV ±1 kV differential mode ±2 kV common mode	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % U ; 0, 5 cycle At 0 , 45 , 90 , 135 , 180 , 225 , 270 and 315 0 % U ; 1 cycle and 70 % U ; 25/30 cycle Single phase: at 0°	N/A	Mains power quality should be that of a typical commercial or hospital environment. If the user of the EQUIPMENT or SYSTEM requires continued operation during power mains interruptions, it is recommended that the EQUIPMENT or SYSTEM be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	N/A	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note : At 80 MHz, the higher frequency range applies.

Note : These guidelines may not apply in all structures. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

- a) Field strength from fixed transmitters, such as base stations for radio (cellular /cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in which the **Oxygo® (DB12)** is used exceeds the applicable RF compliance level above, the monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the monitor.
- b) Over the frequency range 150 kHz to 80 MHz, field strength should be less than [V1] V/m.
- c) The clinical test subject is 12 healthy adult volunteer subjects include 6 men and 6 women (ages 21-41 years, 47-82 KG, height 155-185 cm : 3 volunteers are african, 3 volunteers are caucasian and the others are asian.) With skin color light to dark pigmentation.

Chapter 9 Symbol Definition

Symbol Definition Information

Symbol	DEFINITION	Symbol	DEFINITION
	MANUFACTURER		BF type : Body floating type, higher degree of protection against electronic shock than TYPE B applied part.
	AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY		WEEE SYMBOL, ELECTRONICS TO BE COLLECTED AND DISPOSED SEPARATELY (2012/19/EU)
	PLEASE READ THE USER MANUAL FIRST		ONE LR03 BATTERY, 1.5V DC
	CE MARKING: CONFORM TO EC DIRECTIVE NO.93/42/EEC FOR MEDICAL DEVICES		POWER SWITCH
	CAUTION		NON-CONTINUOUS MONITORING (NO ALARM FOR SPO2)
	SERIAL NUMBER		Solid particle protection: >12.5mm (effective against) Liquid ingress protection: Dripping water when tilted at 15° (effective against)
	STORAGE TEMPERATURE HIGH LIMIT 70°C STORAGE TEMPERATURE LOW LIMIT -25°C		STORAGE RELATIVE HUMIDITY 10% ~ 95% NON-CONDENSING
	ATMOSPHERIC PRESSURE LOW LIMIT 50 KPA ATMOSPHERIC PRESSURE HIGH LIMIT 106 KPA		TO INDICATE MALFUNCTION OR TO IDENTIFY THE CONTROL BY WHICH THE OPERATOR CAN ASK FOR ASSISTANCE[HELP BUTTON]

The product comply with the requirements of the EU directive 2011/65/EU relating to the restrictions of hazardous substances (ROHS). REACH (Registration, Evaluation, Authorization and Restriction of Chemicals, EC 1907/2006) is the European Union's (EU) chemical substances regulatory framework, REACH requires manufacturer to provide sufficient information on Substances of Very High Concern (SVHC).

This product do not contain substances listed on the SVHC which in a concentration above 0.1% weight by weight.

Chapter 10 Customer Service Information

Warning: **Oxygo® (DB12)** Pulse Oximeter is a delicate device and call only be repaired by professional technicians. Disassembling the device without professional training may lead to damage of it.

Within two years from the purchase of this device, SPENGLER will be responsible for the repair should any defects in the material or manufacturing occur. Within two years of the warranty, should the **Oxygo® (DB12)** Pulse Oximeter fail to function properly due to faulty materials or manufacturing, SPENGLER shall replace it with a brand new unit or offer to exchange it with another product of the same value.

? However, failures resulting from personal causes are excluded, such as:

- a) Failure due to modifications of the internal components or disassembling the product without authorization
- b) Failure arising from being dropped during use or in transit
- c) Failure due to the lack of proper maintenance
- d) Failure due to not following the instructions of the User Manual for accurate operations

For any dissatisfaction with the measurement outcomes of the **Oxygo® (DB12)** Pulse Oximeter, please contact the local authorized dealer for more information.

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Imported and distributed by SPENGLER SAS
Delta Electronics (Netherlands) B.V.