Spengler

neotens® nomad

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CONTENTS

Intended Use	47
Contraindications	47
Side Effects	47
Safety Notice	47
Important Instructions before Use	50
Monitor Unit	52
Important Testing Guidelines	54
Quick Start	55
Unit Operation	57
Battery Installation	57
Start Up	58
System Settings	59
Applying the Wrist Monitor	61
Memory Group Selection	62
Testing	62
Irregular Heartbeat Indicator	64
Power Off	64
Memory Check	65
Memory deletion	66
Low Battery Indicator	66
Troubleshooting	67
Error Description	68
Blood Pressure Information	70
Blood Pressure Q&A	72
Maintenance	73
Specifications	76
Materiovigilance	79
Disposal	79
Warranty	79
Mobile Applications	80
Electromagnetic compatibility information.....	81

Thank you for purchasing the Neotens® Nomad (NEO100) Blood Pressure Monitor. The unit has been constructed using reliable circuitry and durable materials. Used properly, this unit will provide years of satisfactory use.

INTENDED USE

This device is intended for non-invasive measurement of systolic and diastolic blood pressure and heart rate in adults and adolescents over 12 years of age using the oscillometric method. The device is not intended for use on infants and children. The device is designed for home or clinical use. All functions can be used safely and values can also be read out on the LCD DISPLAY. The measurement position is placed only on the wrist of an adult. This device is intended for patients. The intended users (also called operators) are also doctors, nurses or care-givers.

CONTRAINDICATIONS

- The cuff should not be applied over damaged skin.

SIDE EFFECTS

- Individuals with serious circulation problems may experience discomfort.
- Too frequent measurements can cause injury to the patient due to blood flow interference.
- Prolonged over-inflation of the cuff may cause bruising of the wrist.
- DO NOT attach the wrist cuff to a limb being used for IV infusions or any other intravascular access, therapy or an arteriovenous (AV) shunt. The cuff inflation can temporarily block blood flow, potentially causing harm to the patient.

SAFETY NOTICE

Please read this manual thoroughly before using the unit. Please retain this manual for future reference. For specific information about your blood pressure, please CONSULT YOUR DOCTOR.

To avoid risk and damage follow all warning precautions. Operate unit only as intended.

WARNING SIGNS AND SYMBOLS USED

	Caution		Keep dry
	Mandatory		Keep away from sunlight
	Prohibited		Power source: batteries
	Type BF equipment		Direct current
	User manual, MUST be consulted		Medical Device
	Serial Number		Do not use if package is damaged
	Batch code		Protected against solid bodies greater than 12.5 mm. Protected against vertical water drops when the shade is tilted up to 15°
	Discard the used product to the recycling collection point according to local regulations.		Temperature limit
	The product conforms to the requirements of the Regulation (EU) 2017/745 MDR on medical devices.		Humidity limitation
	Manufacturer		Unique Device Identifier
	Date of manufacture		Catalogue number
	Authorized representative in Switzerland		Importer
	Authorized representative in the European Community		Distributor
	Model number		Triman for recyclable packaging in France

 Caution

Contact your doctor if test results regularly indicate abnormal readings. Do not attempt to self-treat these symptoms without consulting your doctor first.

Product is designed for its intended use only. Do not misuse in any way.

Product is not intended for infants or individuals who cannot express their intentions.

Do not disassemble or attempt to repair.

Do not use devices, which generate strong electrical or electromagnetic fields near the device, as they may cause incorrect readings and interference or become interference source to the device.

 Battery Precautions:

Do not mix new and old batteries simultaneously.

Replace batteries when Low Battery Indicator "" appears on screen. Replace all batteries at the same time.

Do not insert the batteries with their polarities incorrectly aligned.

Do not mix battery types. Long-life alkaline batteries are recommended.

Remove batteries from device when not in operation for more than 3 months.

Dispose batteries properly; observe local laws and regulations.


IMPORTANT INSTRUCTIONS BEFORE USE

1. Do not confuse self-monitoring with self-diagnosis. Blood pressure measurements should only be interpreted by a health professional who is familiar with your medical history.
2. Contact your doctor if test results regularly indicate abnormal readings.
3. If you are taking medication, consult your doctor to determine the most appropriate time to measure your blood pressure. NEVER change a prescribed medication without first consulting your doctor.
4. For persons with irregular or unstable circulation resulting from diabetes, liver disease, arteriosclerosis or other medical conditions, there may be variations in blood pressure values measured at the wrist versus at the upper arm. Monitoring the trends in your blood pressure taken at either the arm or the wrist is nevertheless useful and important.
5. People suffering from vascular constriction, liver disorders or diabetes, people with cardiac pacemakers or a weak pulse, and women who are pregnant should consult their doctor before measuring their blood pressure themselves. Different values may be obtained due to their condition.
6. People suffering from arrhythmias such as atrial or ventricular premature beats or atrial fibrillation only use this blood pressure monitor in consultation with your doctor. In certain cases, oscillometric measurement method can produce incorrect readings.
7. The wrist cuff should not be placed on the arm on the side of a mastectomy. In the case of a double mastectomy use the side of the least dominant arm.
8. Pressurization of the cuff can temporarily cause loss of function of simultaneously used monitoring equipment on the same limb.
9. Check that operation of the unit does not result in prolonged impairment of the circulation of the patient.
10. The system might produce incorrect readings if stored or used outside the manufacturer's specified temperature and humidity ranges. Make sure to keep the blood pressure monitor away from children or pets.
11. Use only the approved wrist cuff for this unit. Use of other wrist cuffs may result in incorrect measurement results.
12.  Advising operator that the user manual must be consulted.
13. Do not use the device during transport vehicles for influencing measurement accuracy, such as patient transport in an ambulance or helicopter.


14. Contains small parts that may cause a choking hazard if swallowed by infants.
15. The ME EQUIPMENT shall not be used in operating rooms, nor under emergency circumstances, nor during patient transportation or movement.
16. The ME EQUIPMENT shall be usable in general-purpose offices such as professional healthcare facilities.
17. The ME EQUIPMENT shall be usable at home for general-purpose exams.
18. The 5V DC USB connector must not be used. This connector is for qualified service operator only.
19. The device should not be used or exposed to conditions outside its intended specifications.
20. The operator should maintain basic safety and essential performance of the product in accordance with IEC standards.

MONITOR UNIT

Blood Pressure Monitor

Display

Contents

1. Wrist Monitor Unit

2. User manual

3. Batteries AAA (x2)

IMPORTANT TESTING GUIDELINES

1. Avoid eating, exercising, and bathing for 30 minutes prior to testing.
2. Sit in a calm environment for at least 5 minutes prior to testing.
3. Do not stand while testing. Maintain a comfortable sitting position, back and arm supported, legs uncrossed, feet flat on the ground while holding your arm at the level of the heart.
4. Avoid speaking or moving body parts while testing.
5. While testing, avoid strong electromagnetic interference such as microwave ovens.
6. Wait for 1 minute or longer before re-testing.
7. Try to measure your blood pressure at the same time each day for consistency.
8. Test comparisons should only be made when monitor is used on the same arm, in the same position, and at the same time of day.
9. This blood pressure monitor is not recommended for people with severe arrhythmia.
10. Do not use this blood pressure monitor if the device is damaged.

QUICK START

EN

1. Install batteries. (See Figure A)


Figure A

2. Remove clothing from the wrist area. (See Figure B)
3. Rest for several minutes prior to testing. Wrap cuff around left wrist. (See Figure C)

Figure B

Figure C

4. Sit in a quiet place, firmly resting on your chair, with your back supported, legs uncrossed, your feet flat on the floor with the cuff at heart level. (See Figure D)

5. Press "  " button for at least 2 seconds and release it to start the unit. Press the "  " START button again to start the measurement. (See Figure E)

Figure D


Figure E


Note:

The applied part of the ME equipment consists of the cuff which is in direct contact with the patient's skin for a temporary use i.e. continuous usage up to less than 60 minutes.

Battery Installation

Slide battery cover off as indicated by arrow.

Install 2 new AAA alkaline batteries according to polarity. Close battery cover.

Note:

- 1) Replace batteries when "B" blinks on screen.
- 2) Batteries should be removed from device when not in operation for an extended period of time.
- 3) If the batteries inserted in the device appear damaged or leak, the user shall take precautions before touching them.

Start up

1. Power On


Press and hold "  " button for more than 2 seconds then release the button to turn on the device. A beep ^{START} is emitted. The LCD screen is displayed for 1 second when the device performs a quick diagnosis.

NOTE: The user should systematically check that all LCD segments light up and that one short beep is emitted when starting the device. In the event of a defect, the product shall be returned for immediate service.


2. Display Revision

The LCD screen continues to display the hardware and firmware revision identifiers on the SYS and DIA segments, respectively, for one second. For example, the following image shows that hardware and firmware revisions are V1.0 and V1.0.

3. Steady state

The unit continues to enter into steady state, and it wait for user pressing any key to start a new operation.


System Settings

When unit is displaying the basic screen or the measurement results, press  to activate System Settings.

1. Time/Date setting:

"Yr" button and the year figures below (for example: 23 for 2023). To set the year, use the  symbol for (-) and the  symbol for (+).

To validate and change to month, use the  symbol, renew the setting of the month and day. Validate with the  symbol. To set the time system to 12h or 24h, select the one of your choice using the  or  icons. Validate with the  symbol and set the time as you wish in the same way.

2. Save Settings:


While in any setting mode, press "  " button to turn the unit off. All information will be saved.

Note: If unit is left on and not in use for 1 minute, it will automatically save all information and shut down.

Applying the Wrist Monitor

Do not apply over clothing. If wearing a long sleeved shirt, be sure to roll sleeve back to forearm.

Apply the blood pressure monitor to wrist as illustrated. Tighten cuff firmly as not to wiggle.

Note: The blood pressure monitor shall be removed from the wrist in case of patient discomfort due to compression.


Do not stand while measuring. Sit in a comfortable position with your back supported and your feet flat on the floor, with legs uncrossed. Position the center of the cuff at the level of the right atrium of the heart, with the arm supported.

UNIT OPERATION

Memory Group Selection

When the unit is in steady state or displaying a measurement result, you can store the results in two different patient groups by pressing the "P" button, allowing you to switch between the two groups. Test results will automatically store in each selected group.


Testing

When unit is in steady state or displaying measurement result, short press  button for less than 1s to launch measurement(s) according to the selected Patient. When entering the measurement state, the unit will emit 1 beep.

1. Testing


During cuff inflation, the air will slowly increase as indicated by the cuff pressure value. When the pressure is greater than 15 mmHg, the cuff icon  will be displayed.


Note: Keep relaxed while testing. Avoid speaking or moving body parts.

2. Result Display

When normally exiting the measurement state, the unit will emit two consecutive beeps. The screen will display measurements for systolic and diastolic blood pressure, pulse value. When the value of PP (systolic pressure - diastolic pressure) is greater than or equal to 60 mmHg, the PP value will be displayed on the LCD, the mmHg segment will be lit, and the PP icon will blink.

No PP value displayed

PP value displayed

Irregular Heartbeat Indicator

If the blood pressure monitor detects an irregular heart rhythm two or more times during the measurement process, the irregular heart rhythm symbol  appears on the screen with the measurement results. An irregular heart rhythm is defined as a rhythm that is 25% slower or faster than the average rhythm detected during measurement of systolic and diastolic blood pressures. Consult your doctor if the Irregular Heartbeat Symbol  frequently appears with your test results.

Power Off


The " START" button can be pressed for more than 1 second to turn off the device in any mode. The device may turn itself off after approximately 1 minute of inactivity in any mode.

Safety Precaution: If pressure in the wrist cuff becomes too extreme while testing, press the any one of , "", "", ", "START", or "" buttons to stop inflation. off the unit. Beeping sounds are emitted. The cuff pressure will rapidly dissipate once the unit is off.

Memory Check



When the unit is in steady state, you may check past test results by using the " " button. Upon activating, test results can be scrolled through using the " " and " " buttons to browse stored results. The LCD below will display the last measurement stored in memory: e.g. measurement 01.

Note: In the event of a measurement data shown as "EEE", the product shall be return for maintenance immediately.

Memory deletion

All memory may be deleted while in Memory Check mode. Press and hold "⏏" for about 3 seconds to delete all memory records for all patient groups. The LCD screen displays "- - -" and beeps.


Note: Memory cannot be recovered once it has been deleted.

Low Battery Indicator

The unit will display E08 when battery life is depleting and unable to inflate cuff for testing.

When "⏏" flashes, it means the battery is lower than 25% and the cuff cannot be inflated for measurement.

The "⏏" appears simultaneously for approximately 5 seconds prior to shutting off. Replace batteries at this time. No memory loss will occur throughout this process.

Troubleshooting

EN

Anomaly	Possible cause	Solution
Blood pressure monitor working abnormally	The wrist cuff is too tight or too loose, or the strap wrist is poorly attached	Roll the cuff correctly
	Arm/wrist or device movement during measurement of blood pressure	Stay quiet, keep your arm/wrist steady and do not move the monitor.
	Talking, being nervous or emotional during measurement	Avoid speaking, breathe deeply to calm down and try to relax.
	Incorrect measurement posture	Adjust your posture.
	There is interference in the loading process or an incorrect operation in the measurement process.	See operating instructions.

The following table shows the error messages that may occur during the measurement, possible causes and treatment methods. Please measure again using the correct method.

When an error or warning occurs, the error code is displayed on the screen and the cuff icon flashes.

The type of code is a warning: one beep is emitted, short press the "  " button to return to the stable state.

The type of code is an error: emits 4 continuous beeps. Displays the error 3 seconds, then automatically powers off.

Error Description

Code	Type	Description of Issue	Possible Causes and Solutions
E01	Warning	The device is not working properly.	Place the device on the wrist again and check that it is working properly.
E02	Error	The deflation time of the cuff is too long.	Please place the device on the wrist again, perform a new measurement and check if the pressure is correct
E03	Error	The pressure detected exceeds 300 mmHg	Please place the device on the wrist again, perform a new measurement and check if the pressure is correct.
E04	Warning	Abnormal measurement data	Please check the wrist cuff is worn properly before re-measuring.
E05	Warning	The inflation speed of the wrist cuff does not meet standards	Please check the wrist cuff is worn properly before re-measuring.
E06	Warning	Cannot complete the measurement normally	Please check the wrist cuff is worn properly before re-measuring.
E07	Warning	Measurement data exceed the standard range	Please readjust the wrist cuff and take another measurement. If the problem is not solved, please contact the manufacturer.
E08	Error	Low battery	Replace batteries.

E09	Warning	The pressure exceeds 15 mmHg for more than 3 minutes	After the warning , the arm will automatically deflate and a new measurement can be taken.
E10	Warning	Arm shaking when measuring	Hold the arm or body still and measure again.
E20	Error	Firmware check error	Please contact the manufacturer.
E21	Error	Abnormal internal voltage	Please contact the manufacturer.
E41	Warning	Pressure offset check error.	Please contact the manufacturer.


Note: If you cannot resolve the abnormal situation by yourself, you can consult Spengler Technical support center or your retailer. It is forbidden to disassemble and repair without permission.

BLOOD PRESSURE INFORMATION

Blood pressure is the force of blood pushing against the walls of arteries. It is typically measured in millimeters of mercury (mmHg). Systolic blood pressure is the maximum force exerted against blood vessel walls each time the heart beats. Diastolic blood pressure is the force exerted on blood vessels when the heart is resting between beats.

An individual's blood pressure frequently changes throughout the course of a day. Excitement and tension can cause blood pressure to rise, while drinking alcohol and bathing can lower blood pressure. Certain hormones like adrenaline (which your body releases under stress) can cause blood vessels to constrict, leading to a rise in blood pressure.

If these measuring numbers become too high, it means the heart is working harder than it should.

Example: fluctuation within a day (male, 35 years old)

BLOOD PRESSURE INFORMATION

The European Society of Cardiology (ESC) defined different levels of blood pressure based on clinical trials and meta-analyses. The prerequisite of Blood Pressure measurement is that it must be undertaken using an accurate device that has been clinically validated. The blood pressure measurement and its monitoring can be performed at patient home; however, the hypertension diagnoses will be established by the doctor at his/her office. According to the environment where the measurement is made, the thresholds test values for elevated pressure or hypertension may vary. The ESC proposed table below summarizes those thresholds at the doctor office or at home.

Non-elevated blood pressure	Elevated blood pressure	Hypertension
Office BP SBP < 120 mmHg and DBP < 70 mmHg	Office BP SBP 120 - 139 mmHg or DBP 70 - 89 mmHg	Office BP SBP ≥ 140 mmHg or DBP ≥ 90 mmHg
Home BP SBP < 120 mmHg and DBP < 70 mmHg	Home BP SBP 120 - 134 mmHg or DBP 70 - 84 mmHg	Home BP SBP ≥ 135 mmHg or DBP ≥ 85 mmHg
Insufficient evidence confirming the efficacy and safety of BP pharmacological treatment.	Risk stratify to identify individuals with high cardiovascular risk for PB pharmacological treatment.	Cardiovascular risk is sufficiently high to merit BP pharmacological treatment initiation.

Source: Blood pressure categories. McEvoy, et al. 2024 ESC Guidelines for the management of elevated blood pressure and hypertension. Eur. Heart Journal.

BP, blood pressure; DBP, diastolic blood pressure; SBP, systolic blood pressure.

Note: Do not be alarmed if an abnormal reading occurs. A better indication of an individual's blood pressure occurs after 2-3 readings taken at the same time each day over an extended period of time. Consult your doctor if test results remain abnormal.

BLOOD PRESSURE Q&A

Q: What is the difference between measuring blood pressure at home or at a professional healthcare clinic?

A: Blood pressure readings taken at home are now seen to give a more accurate account as they better reflect your daily life. Readings can be elevated when taken in a clinical or medical environment. This is known as White Coat Hypertension and may be caused by feeling anxious or nervous.

Note: Abnormal test results may be caused by:

1. Improper cuff placement

Make sure cuff is snug, not too tight or too loose.

Make sure bottom of the cuff is approximately 1-2 cm (1/2") away from the wrist joint.

2. Improper body position

Make sure to keep your body in an upright position.

3. Feeling anxious or nervous

Take 2-3 deep breaths, wait a few minutes and resume testing.

Q: What causes different results?

A: Blood pressure varies throughout the course of a day. Many factors including diet, stress, cuff placement, etc. may affect an individual's blood pressure.

Q: Should I apply the cuff to the left or right arm? What is the difference?

A: Either arm can be used when testing, however, when comparing results, the same arm should be used. Testing on your left wrist may provide more accurate results as it is located closer to your heart.

Q: What is the best time of day for testing?

A: Morning time or any time you feel relaxed and stress free.

MAINTENANCE


EN

1. Avoid dropping, slamming, or throwing the unit.

The user must ensure that the device is positioned on a stable surface. In case the device drops, it may have sharp parts. In case of breakage, the operator must discontinue using the product and follow the instructions outlined in the "Disposal" section.


2. Avoid extreme temperatures. Do not expose unit directly under sunshine.

3. When cleaning the unit, use a soft fabric and lightly wipe with a mild detergent. Use a damp cloth to remove dirt and excess detergent. Disinfection products may include WIP'ANIOS EXCEL, SURFA'SAFE Premium, ANIOSPRAY QUICK, ANIOSPRAY SURF 29. Use a soft cloth to wipe. Only use recommended cleaning and disinfecting products and procedures.

4. Cuff Cleaning: Do not soak cuff in water! Apply a small amount of rubbing alcohol to a soft cloth to clean cuffs surface. Use a damp cloth (water based) to wipe clean. Allow cuff to dry naturally at room temperature. The cuff must be cleaned and disinfected before use between different users.

5. Do not use gasoline, diluents or similar solvents.

6. Remove batteries when not in operation for an extended period of time.

EN

7. Do not disassemble product.

8. It is recommended that the performance is checked every 2 years.

9. Expected service life: Approximately three years at 10 tests per day.

10. Service and maintenance must not be performed while the device is in use. The maintenance can only be performed by service personnel. Service and maintenance require parts, repair and technical support.

Contact Spengler technical support center by phone: +33 (0)4 42 90 31 31 or by e-mail: SAV@spengler-med.fr

SPECIFICATIONS

Product Description	Automatic wrist blood pressure monitor	
Model	NEO100	
Display	LCD Digital Display Size: 32.5 x 40.5mm (1.28" x 1.59")	
Measurement Method	Oscillometric method with measurement during inflation	
Measurement Range	Systolic Pressure	60 mmHg-260 mmHg
	Diastolic Pressure	40 mmHg-200 mmHg
	Pressure	0 mmHg-299 mmHg
	Pressure accuracy	±3 mmHg
	Pulse	30 ~ 199 Beats/Minute
	Pulse accuracy	±5%
Pressurization	Automatic Pressurization	
Memory	50 Memories per user (x2) with Date and Time	
Function	Irregular Heartbeat Detection	
	Low Battery Detection	
	Beep	
	Smartphone application compatibility	
	Automatic Power-Off	
Power Source	2 Alkaline Batteries Size AAA 1.5V	
Battery Life	Approx. 300 measurements with new batteries	
Unit Weight	Approx. 96g (Excluding Batteries)	
Unit Dimensions	Approx. 70.9 ×70.9 ×26.5 mm(L x W x H) (2.79" x 2.79" x 1.04")	

Cuff Circumference	Fits wrist circumference 13.5-21.5 cm(5.3"-8.5")		EN
Operating Environment	Temperature	5°C ~ 40°C (41°F-104°F)	
	Humidity	15%-93%RH	
	Pressure	700 hPa - 1060 hPa	
Storage and Transport Environment	Temperature	-25°C ~70°C (-13°F-158°F)	
	Humidity	≤93% RH	
Ingress Protection Rating	IP 22		
Classification	Type BF equipment  Cuff is the Applied Part		

Specifications are subject to change without notice.

This blood pressure monitor complies with the European regulations and bears the CE mark "CE 0123".

This blood pressure monitor also complies with mainly following standards (included but not limited):

Safety standard:

EN 60601-1 Medical electrical equipment part 1: General requirements for basic safety and essential performance

EN 60601-1-2 Medical Electrical Equipment -- Part 1-2: General Requirements For Basic Safety And Essential Performance -- Collateral Standard: Electromagnetic Disturbances – Requirements And Tests

Performance standards:

IEC80601-2-30, Medical electrical equipment – Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers.

ISO 81060-2, non-invasive sphygmomanometers - part 2: clinical validation of automated measurement type.

MATERIOVIGILANCE

Any serious incident occurring in connection with the device must be notified to the manufacturer and the competent authority of the member state where the user is established.

DISPOSAL

(Waste Electrical & Electronic Equipment)

This marking shown on the product indicates that it should not be disposed of with other household waste at the end of its life. To prevent potential harm to the environment or to human health, please separate this product from other types of wastes and recycle it responsibly. When disposing this type of product, contact the retailer where product was purchased or contact your local government office for details regarding how this item can be disposed in an environmentally safe recycling center.

Business users should contact their supplier and check the terms and conditions of the purchasing agreement. This product should not be mixed with other commercial wastes for disposal. This product is free of hazardous materials.

WARRANTY

The blood pressure monitor is guaranteed for 3 years from the date of purchase. If the blood pressure monitor does not function properly due to defective components or poor workmanship, we will repair or replace it freely. In this case, contact technical support center by phone 04 42 90 31 31 or email sav@spengler-med.fr. The warranty does not cover damages to your blood pressure monitor due to improper handling. Please contact local retailer for details.

MOBILE APPLICATIONS

The blood pressure monitor is compatible with the **Spengler HTA** application which allows you to take a photo of the results or to enter them directly.

The blood pressure monitor is also compatible with the **SuiviHTA** application of the FRHTA. Scan the code below to access the applications.

ELECTROMAGNETIC COMPATIBILITY INFORMATION

EN

The device satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

Table 1

Guidance and declaration of manufacturer-electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment -guidance
Radiated emission CISPR 11	Group 1, class B.	The device uses RF energy only for its internal function. Therefore, its emissions are very low and are not likely to cause any interference in nearby electronic equipment
Conducted emission CISPR 11	N/A	
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	N/A	

Table 2

Guidance and declaration of manufacturer-electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment -guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient/ burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	N/A	
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	N/A	

Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycle 70% UT (30% dip in UT) for 25 cycle <5% UT (>95% dip in UT) for 5 cycle	N/A	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m ; 50 Hz or 60Hz	30 A/m ; 50 Hz or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Radiated RF EM fields IEC 61000-4-3	3V/m or 10 V/m 80MHz-2.7 Ghz 80%AM at 1kHz	3V/m or 10 V/m 80MHz-2.7 Ghz 80%AM at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance 80 MHz to 800 MHz 800 MHz to 2.7 Ghz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an Electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Conducted disturbances Induced by RF fields IEC 61000-4-6	3 V in 0.15 MHz- 80 MHz 6 V in ISM and/or amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1kHz	3 V in 0.15 MHz- 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1kHz	

Table 3

Guidance and declaration of manufacturer-electromagnetic immunity						
Nowadays, many RF wireless equipments have being used in various healthcare locations where medical equipment and/or systems are used. When they are used in close proximity to medical equipment and/or systems, the medical equipment and/or systems' basic safety and essential performance may be affected. Arm-type Fully Automatic Digital Blood Pressure Monitor has been tested with the immunity test level in the below table and meet the related requirements of IEC 60601-1-2:2014+AMD1:2020. The customer and/or user should help keep a minimum distance between RF wireless communications equipment and this medical equipment and/or systems as recommended below.						
Test frequency (MHz)	Band (MHz) ^{a)}	Service ^{a)}	Modulation	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/M)
385	380-390	TETRA 400	Pulse Modulation 18Hz ^{b)}	1.8	0.3	27
450	430-470	GMRS 460 FRS 460	FM $\pm$ deviation de 5 kHz 1 kHz sine	2	0.3	28
710	704-787	LTE Band 13, 17	Pulse Modulation 217Hz ^{b)}	0.2	0.3	9
745						
780						

810	800-960	GSM 800 /900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse Modulation 18Hz ^{b)}	2	0.3	28
870						
930						
1720	700-1990	GSM 1800; CDMA 1900 ; GSM 1900 ; DECT ; LTE Band 1, 3, 4, 25; UMTS	Pulse Modulation 217Hz ^{b)}	2	0.3	28
1945						
1970						
2450	2400-2570	Blue tooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse Modulation 217Hz ^{b)}	2	0.3	28

5240	5100-5800	WLAN 802.11 a/n	Pulse Modu- lation 217Hz ^{b)}	0.2	0.3	9
5500						
5785						
If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.						
a) For some services, only the uplink frequencies are included. b) The carrier shall be modulated using a 50 % duty cycle square wave signal. c) As an alternative to FM modulation, the carrier may be pulse modulated using a 50 % duty cycle square wave signal at 18 Hz. While it does not represent actual modulation, it would be worst case.						

Table 4

Recommended separation distances between portable and mobile RF communications equipment and the device		
The device is intended for use in an electromagnetic environment in which radiated therefore disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.		
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)	
	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.7 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.23
0.1	0.38	0.73
1	1.2	2.3
10	3.8	7.3
100	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

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